

INDATECH READY

INDATECH's multitasking software



Control of all Indatech instruments

Importation of chemometric models

Simplified mode for operators


Communication via OPC-UA

Validated as CFR21 part 11

For the biotechnological and pharmaceutical industries, Indatech - Chauvin Arnoux Group proposes several innovative spectroscopic analyzers for in-line assessment of your product's critical physical and chemical quality criteria in real time. From upstream to downstream, and not forgetting the lyophilization stage, our experts offer you the best possible cost/performance trade-off by taking your objectives into account: speeding up development of the process or improving and automating production. Our solutions, which are fully compatible with automation thanks to their OPC UA communication capabilities, have already been integrated in the systems of the major players such as GEA, Pharmatechnology and Proditec.

Multitasking software.

Software for checking, collecting and processing instantaneous data transmitted by the measuring instrument, from the laboratory through to the production line.



Adapted for each INDATECH solution

Software versions for the whole range of INDATECH spectroscopy instruments:

Asuryan Ready (UV VIS), **Viserion Ready** (Raman), **SRS Ready** (NIR with SAM SPEC technology).

Data processing

Ability to use the full power of chemometric analysis:

- Possibility of directly importing models from Eigenvector Research's PLS-Toolbox,
- Ability to work with the SIMCA suite from Sartorius.
- Display of the predictions in "READY" in real time.

State at delivery

The Ready software is delivered with all the instruments in the Indatech range to allow control of the instrument and interaction with the developers, operators or other machines.

Control by automation

The software includes advanced communication protocols such as OPC-UA or TCP IP. In this way, it is possible to control and command the instruments remotely by means of automation or via advanced production software such as SIPAT, SyntQ or SCADA. Indatech's instruments thus become totally autonomous on a production unit, saving time for the operator.



21 CFR Part 11

Compatible with pharma standards

Acquisition in an internal data format (IDM "Indatech Data Measurement"), or the more standard "SPC" format, both suitable for spectroscopic measurement, with speed varying according to the technology used. Tools are also provided for instantaneous extraction of the data and activity reports as well as **management of the access rights or electronic signatures**, all in conditions compliant with the **21 CFR part 11 standard**.

User	Access Level	Status	Creation Date	Password Update	Max Login Attempts	Password Lifetime (MONTHS)
OM	INDATECH	OK	06/11/2020	20/11/2020	10	24
YT	INDATECH	OK	06/11/2020	16/11/2020	10	24
FC	INDATECH	Expired	06/11/2020	01/01/1904	10	24
SBA	INDATECH	Blocked	06/11/2020	16/11/2020	10	24
LA	INDATECH	OK	06/11/2020	23/11/2020	10	24
DB	INDATECH	OK	06/11/2020	10/06/2021	10	24
JL	INDATECH	OK	06/11/2020	26/11/2020	10	24
CG	INDATECH	OK	01/01/1904	03/06/2022	10	24
TS	INDATECH	OK	01/01/1904	26/11/2021	10	24
RP	INDATECH	Expired	01/01/1904	01/01/1904	10	24
AdminUser	Administrator	OK	27/12/2021	03/06/2022	3	12
SecondAdmin	Administrator	OK	27/12/2021	27/12/2021	3	12
AutoTest	User	Expired	09/06/2022	09/06/2022	1	7

Data management

- A system for modulating the data acquisition parameters:
 - > type of acquisition (one-off, periodic or configured via recipe),
 - > integration time,
 - > number of scans.
- Acquisition in an internal file format containing the various data of interest (references, spectra, calculated process values for the object, metadata, experimental environment, batch no., etc.), or in standard SPC format.
- Possibility of exporting in the universal SPC format
- A data acquisition and recording procedure designed for compatibility with the 21 CFR part 11 standard
- Generation of end-of-batch reports



Generation of recipes for operators

- Definition of the acquisition and display parameters
- Importation of chemometric models (PLS, ACP, SVM, ANN, etc.)
- Definition of the acceptance limits for each of the imported models
- Automatic loading of the parameters for the operator

Human/machine interface

The screenshot displays the INDATECH READY software interface. It features a main window with a title bar and a menu bar. The main area is divided into several sections:

- Title:** The top bar of the application window.
- Menu:** The top-left corner containing options like File, Display, Measurement, Settings, Setup, Tools, and System Validation.
- System status:** The top-right corner showing system information and a Log button.
- Configuration panel:** A vertical panel on the right side with settings for three channels (Ch 1, Ch 2, Ch 3), including integration time, scans, and buttons for 'Set as Reference', 'Save References', and 'Save Last Spectra'.
- Activation of acquisition mode:** A 'Process' button at the bottom right of the configuration panel.
- Data display:** The central area containing three plots: 'Spectrum' (Intensity vs Wavelength), 'Spectrum Desk' (Intensity vs Wavelength), and 'Absorbance' (Absorbance vs Wavelength).
- Data backup information:** The bottom status bar showing the current directory and file name.

Data display

- Possibility of viewing and monitoring the intensity, dark and absorbance in real time.
- Possibility of displaying a history of the most recent spectra
- Possibility of multiscreen display
- Display of the critical parameters for the chemometric models (prediction, Hotelling's, Q^2).


Qualification, maintenance and control tools

- Qualification of the equipment according to the applicable USP standards
- A tester of the status of the various elements with a limited life span (light source)
- A tool for routine verification of the system with an alarm for checking the qualification of the equipment (spectrophotometers, light source, etc.), and your experiments (dark, reference)
- Export of "daily check" reports

To order

Measuring system	Type of technology	USP qualification	21 CFR part 11	Importation of models		Communication mode		
				PLS Toolbox	Simca (option)	OPC-UA (option)	TCP-IP	I/O
Asuryan	UV-Visible	USP 857	Compatible	Yes	Yes	Yes	No	No
Hyternity	NIR	USP 856	Compatible	Yes	Yes (low speed)	Yes (low speed)	Yes	Yes
Viserion	Raman	USP 858	Compatible	Yes	Yes	Yes	No	No

Multipurpose software



906211679 - Ed. 1 - LCC - 12/2022 - Non-contractual document.

CA INDATECH

4, rue Georges BESSE - 34 830 Clapiers - France
 Tel.: +33 4 80 78 01 40 | Fax: +33 4 80 78 01 42
 sales@indatech.eu | www.indatech.eu

